SILHOUETTE INSTALIFT™

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SMS 31

Instructions for Use
Directives d'utilisation

Manufactured for:
Sinclair Pharma US Inc.
18 Technology Drive
Unit 134, Irvine,
CA 92618 USA

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PN: SMS 31-PLG-3.0.1-NA

CONTENTS OF THE CARTON:

SMS 28: 5 packs of SMS 28-PLG-3.0.1-NA (2 units), 5 x patient chart labels, 1 x Ruler, 1 x Instructions for Use.

SMS 29: 5 packs of SMS 29-PLG-3.0.1-NA (2 units), 5 x patient chart labels, 1 x Ruler, 1 x Instructions for Use.

SMS 30: 5 packs of SMS 30-PLG-3.0.1-NA (2 units), 5 x patient chart labels, 1 x Ruler, 1 x Instructions for Use.

SMS 31: 5 packs of SMS 31-PLG-3.0.1-NA (2 units), 5 x patient chart labels, 1 x Ruler, 1 x Instructions for Use.

INDICATIONS:

Silhouette Instalift is a resorbable sterile implantable single-use device intended for adult patients and is used in multiple pairs to provide for elevation and repositioning of facial tissue.

DESCRIPTION:

The Silhouette Instalift device is a resorbable implantable suspension suture for facial tissue lifting, consisting of a resorbable sterile monofilament USP 3-0/ EP 2 size and injection molded cones affixed to it. Both components are made of Poly (L-lactide-co-glycolide) 82:18. The 30cm suture (SMS 28-PLG-3.0.1- NA) has 8 cones, 27.5cm suture (SMS 29-PLG-3.0.1-NA/ SMS 31-PLG-3.0.1-NA) has 12 cones and the 26.8cm suture (SMS 30-PLG-3.0.1-NA) has 16 cones. All have two 12cm straight needles attached to each end.

The number of sutures required per treatment is dependent on the degree of repositioning required, the size and anatomy of the treatment area, as well as physician technique and recommendation. A typical treatment will use up to 4 sutures per side.

All products are supplied sterile (EO) for single use only.

ACTIONS:

Silhouette Instalift elicits a minimal acute inflammatory reaction in tissue that is followed by gradual encapsulation. The device is for placement in subcutaneous tissue.

DURATION OF EFFECT:

Silhouette Instalift treatment effects have variable duration* dependant on the insertion pattern used (by the physician) and the patient's lifestyle

**clinical data on file has shown a good duration of effect up to 24 months follow-up*

PRECAUTIONS FOR USE:

The physician should be familiar with recommended techniques involving Silhouette Instalift as well as proper patient selection and device placement. Ensure procedure is carried out in a sterile environment. Knowledge of anatomy of the site to be treated and specific precautions are essential in order to avoid damage to fragile structures such as nerves and blood vessels. Patients with extensive skin laxity, thin skin and/or severe malar fat sagging are not good candidates for this procedure. Take extra precautions for patients with neurological conditions. Take extra care with patients who have bleeding disorders and with patients who are taking medications that inhibit blood clotting (e.g. vitamin E, acetylsalicylic acid and salicylic acid derivatives and nonsteroidal anti-inflammatory drugs).

CONTRAINDICATIONS:

Patients with foreign body sensitivity or known or suspected allergies to implant or instrument materials in particular plastic/biomaterial should not undergo the procedure. Do NOT use in patients with active sepsis or infection, active (or history of) autoimmune disease, patients under 18 years of age, pregnant or breastfeeding women, or patients with limited ability or un-willingness to follow post-treatment recommendations.

INSTRUCTIONS:

1. Ensure patients face is make-up free, and clean the treatment site with antiseptic.
2. Mark the entry and exit locations where the Silhouette Instalift suspension sutures are to be inserted (Fig 1 shows example placement for mid face treatment).
3. Remove sutures from protective pouch and place in sterile field using aseptic procedures.
Warning: Sutures are double pouched and only the internal surface of the inner pouch is sterile.
4. Administer appropriate local anesthetic on entry and exit points (Fig 2 example for mid face treatment).
5. Make a pre-incision hole with an 18G needle (not supplied) in the treatment point marked for insertion.
6. Insert one of the suture needles through the pre-incision hole to the required depth in the subcutaneous tissue as indicated by the black mark on the needle (Fig.3) and advance along the pre-marked location for each suture placement exiting the skin in the desired area (Fig 4). Run the suture through the skin by "pinch traction".
7. The needle exits through the pre-marked exit point and then the needle may be cut off leaving a portion of the suture exposed (Fig 5).
8. Repeat using second needle to insert the other half of the suspension suture (Fig 6-9).
9. Apply tension to the suture in place to achieve the desired lift and volume (Fig 10).
10. Trim the distal ends of the suture. It is desirable to bury the ends under adjacent, mobilized soft tissue to avoid possible erosion of the suture through the skin.
11. Repeat for each Silhouette Instalift suspension suture implanted.

POST-TREATMENT RECOMMENDATIONS:

Apply cold packs immediately after procedure if required (cold packs should be wrapped to avoid direct contact with skin and insertion point). Take analgesia, e.g. paracetamol, in case of pain. Refrain from applying make-up for as long as possible (minimum 24 hours). Sleep face up, in an elevated position for 3-5 nights. Wash, shave and dry the face gently without rubbing or massaging (5 days). Avoid excessive face and neck movements (2 weeks). Avoid over-exposure to direct sunlight and do not use tanning beds (2 weeks). Avoid high impact sports e.g. running (2 weeks). Do not use saunas and hammams (3 weeks). Avoid dental surgery (3 weeks). Avoid face and neck massages and facial aesthetic treatments (4 weeks).

ADVERSE EVENTS:

Like all procedures of this type there is a possibility of adverse events, although not everybody experiences them. These adverse events include, but are not limited to, infection, minimal acute inflammatory tissue reaction, redness, pain (which may be temporary or persistent in nature), tenderness, swelling and oedema (including nodules), transient haematoma or bruising and transient rippling or dimple formation. Other potential adverse events include ecchymosis, sensory/motor nerve injury, asymmetry, banding, puckering, thread migrations, palpable thread ends/knots/cones and visible or protruding threads/cones. Material sensitivity/allergic reactions in patients following surgery may occur. Lack of effect has also been reported. The onset of any adverse event must be reported immediately. Please contact the local Sinclair representative or authorised Silhouette distributor. Alternatively send the details to quality@sinclairpharma.com

WARNINGS:

In the handling of the suspension suture and any other suture-like material, care should be taken to avoid damage from handling, if the device appears damaged prior to use, do not use. Avoid crushing or crimping damage. The superficial placement of this device may lead to all of the aforementioned side effects. Too deep placement of the suture (in the muscle) can cause pain. Do not use if the packaging is damaged or open. Although the combination use of Silhouette devices with other aesthetic procedures has been reported in literature, the safety of such combination treatments has not been formally established.

Single use only. Re-use of this suspension suture may result in hazardous biocontamination resulting in severe injury to the patient. Do not re-sterilize.

Discard open, unused product. Product should be disposed of in a manner that complies with the regulatory requirements of the country where it is sold.

Caution: Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare practitioner.

Store between 0°C and 8°C in a dry place out of direct sunlight.
Do not use beyond the listed expiration date on package.
If you have any complaints then please contact
quality@sinclairpharma.com

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